

## Autoclave deodorisers, capsules 03



Neutralises pungent malodours from autoclaving equipment soiled with organic matter.

Gelatine capsules containing highly concentrated deodorising materials.

Suitable for use in all makes of laboratory, medical and industrial autoclaves.

Catalogue No	Pack qty
AUX-050-010H	100
AUX-050-012D	500

## Autoclave tape, sterilising, indicator 04



Process indicator tape that meets BS 7893: 1997 specification for pressure sensitive adhesive closing and sealing tapes for use with sterilisation packaging materials.

Also meets BS EN 867-21997 non-biological systems for use in sterilisers Part 2. Process indicators (Class A).

Catalogue No	Colour	Width, mm	Pack type
SVF-610-010J	Cream	19	Roll/50m
SVF-610-020G	Blue	19	Roll/50m
SVF-610-030D	Cream	25	Roll/50m

## Autoclave tape labels, indicator 04

Ideal for labelling glassware, can be written on.

Process indicator labels that meet BS 7893: 1997 specification for pressure sensitive adhesive. Do not leave residue when applied to glass, etc.

Also meet BS EN 867-21997 non-biological systems for use in sterilisers Part 2. Process indicators (Class A).

Catalogue No	Description	Width, mm
SVF-610-070X	Roll of 500 labels	19

## Autoclave tape 04



High temperature resistant tape, with special adhesive qualities to withstand steam sterilisation.

Available plain or with steam sensitive indicator ink. The colour change is designed to show that the pack has been exposed to a steam sterilisation process.

Standard tape

Catalogue No	Width, mm	Pack type
AUY-210-010J	19	Roll/50m
AUY-210-020G	25	Roll/50m

Steam sensitive ink indicator tape

Catalogue No	Width, mm	Pack type
AUY-240-010L	19	Roll/50m
AUY-240-020Y	25	Roll/50m

## Autoclave tape, indicator 04


For use with sterilisation packaging materials

Meets BS 7893: 1997 specification for pressure sensitive adhesive closing and sealing tapes, also meets BS EN 867-21997 non-biological systems for use in sterilisers Part 2. Process indicators (Class A).



Catalogue No	Width, mm	Pack type
AUY-170-030P	18	Roll/50m
AUY-170-040M	24	Roll/50m

## Sterilisation indicators, chemical, tubes 02



Effective and immediate visual check on sterilisation conditions.

Integrated 2 criteria control of temperature and time for steam, fluid and dry heat sterilisers.

Catalogue No	Description	Pack qty
AUX-200-010D	Black spot, 121°C/15min <sup>1</sup>	100
AUX-200-020A	Yellow spot, 134°C/3min <sup>1</sup>	100
AUX-200-030U	Green spot, 160°C/60min <sup>1</sup>	100
AUX-200-040H	Blue spot, 180°C/12min <sup>1</sup>	100
AUX-200-050X	White spot, 160°C/120min <sup>1</sup> /170°C/68 min <sup>1</sup> /180°C/35min <sup>1</sup>	100

# Autoclaves

We offer a comprehensive range of portable, benchtop and floorstanding laboratory autoclaves and accessories from all of the leading manufacturers, please contact our sales dept with your specific requirements for further details.

