

Beakers, borosilicate glass, conical form, with spout

01

| Catalogue No | Alt. No | Capacity, mL |
|--------------|----------|--------------|
| BNB-410-110Y | 161/150 | 150 |
| BNB-410-130C | 161/250 | 250 |
| BNB-410-160Q | 161/500 | 500 |
| BNB-410-230V | 161/1000 | 1,000 |



Beakers, Pyrex® glass, conical form, with spout

01

| Catalogue No | Alt. No | Capacity, mL | Diameter [external], mm | Height, mm | Pack qty |
|--------------|----------|--------------|-------------------------|------------|----------|
| BNB-425-030P | 1020/06D | 250 | 68 | 52 | 10 |
| BNB-425-070D | 1020/16D | 500 | 86 | 59 | 10 |



Beakers, polypropylene, squat form, BS 5404 : Part 1, ISO 7056, Thermo Scientific Nalgene

55

Tapered with ribs for easy stacking. Clear graduations. Autoclavable.

| Catalogue No | Alt. No | Capacity, mL | Pack qty |
|--------------|-----------|--------------|----------|
| BNH-680-010C | 1201-0030 | 30 | 12 |
| BNH-680-030T | 1201-0050 | 50 | 12 |
| BNH-680-050N | 1201-0100 | 100 | 12 |
| BNH-680-060K | 1201-0150 | 150 | 12 |
| BNH-680-070H | 1201-0250 | 250 | 6 |
| BNH-680-090B | 1201-0400 | 400 | 6 |
| BNH-680-110V | 1201-0600 | 600 | 4 |
| BNH-680-130P | 1201-1000 | 1,000 | 3 |
| BNH-680-150J | 1201-2000 | 2,000 | 1 |
| BNH-680-170U | 1201-4000 | 4,000 | 1 |

